



# PORTABLE CLEAN ROOM

## SS-330-PCR

Weight: Approx. 70 lbs  
HEPA Filtration - Up to 99.97%

### Features

- Adjustable Speed Control
- Quiet Operation
- Attractive Appearance
- Solid Construction
- Patented Fan/Filter Blower
- Hinged Sloping Lid.
- Low Maintenance Operation
- Long Filter Life

### Typical Applications

- Clean Rooms
- Medical
- Pharmaceutical
- Commercial
- Laboratory
- Sterile Cleaning Environment

### Technical Specifications

#### Hood Dimensions:

O.D: 24.5" D x 32" W x 25" H  
I.D.: 24" D x 30" W x 24" H

#### Hood Material:

Sides: .250" White Acrylic  
Base: .500" White High Density Polyethylene  
Top: .500" Clear Polycarbonate  
Front: .250" Clear Polycarbonate

#### Blower Dimensions:

15" L x 13" W x 13" H

#### Blower Material:

16 ga. carbon steel  
w/ polyester epoxy coating on all surfaces

#### Motor:

Fractional Horsepower, TEFC

#### Electrical:

115/1/60, Approx. 2 amps

#### Sound Level:

Approx. 55 dba @ 3'

When the application calls for a positive pressure laminar flow hood, the Model SS-330-PCR 30"W Portable Clean Room for Positive Pressure may be an excellent choice. The powerful fan on top of the hood combined with a variety of Hepa or Ulpa filter options provides exceptionally clean, purified air-flow inside the hood. Typical uses vary widely, but are predominantly this system is used when an economical and portable benchtop cleanroom environment is called for.

Sentry Air's Portable Clean Room Hoods can provide an effective and economical solution for many commercial, medical, laboratory and industrial applications.

### Other Models Available

- SS-340-PCR - 40"
- SS-330-PCR - 30"
- SS-324-PCR - 24"
- SS-218-PCR - 18"
- SS-212-PCR - 12"

Custom Sizes & Models Available

USA: 800.799.4609  
INTL: 713.690.2153  
FAX: 713.690.7872

**800.799.4609**  
**WWW.SENTRYAIR.COM**

Sentry Air Systems, Inc.  
6999 West Little York RD., #P1  
Houston, TX 77040