


PORTABLE CLEAN ROOM

30" & 40"

MODELS: SS-330-PCR & SS-340-PCR

When the application calls for a **Positive Pressure Flow Hood**, the Sentry Air Systems 30" & 40" Portable Clean Rooms are sure to be an excellent choice. HEPA or ULPA purified air offers superior protection from unwanted particulates.

- » HEPA/ULPA FILTRATION
- » QUIET OPERATION
- » CONVENIENT ACCESS
- » COST EFFECTIVE



The powerful, patented blower, combined with HEPA or ULPA filters, provides purified airflow. Sentry Air Systems' Portable Clean Rooms offer a cost effective and superior alternative in clean air applications.

TYPICAL USES:

- ◆ Clean Rooms
- ◆ Forensic Laboratories
- ◆ Manufacturing Areas
- ◆ Sterile Service Operations
- ◆ Biological Applications

ADDITIONAL FEATURES:

- ◆ Laboratory Safe Materials
- ◆ Durable Construction
- ◆ Hinged Sloping Lid
- ◆ High Airflow
- ◆ Quiet Operation

OPTIONS:

- ◆ Custom Designs Available
- ◆ Hour Counter
- ◆ 27" Fluorescent Light
- ◆ Anti-Static Strip Curtains
- ◆ Clean & Polish Kit

FILTERS: HEPA (Up to 99.97% @ .3 microns), ULPA (Up to 99.9995% @ .12 microns) & Polyester Pre-Filter Pads.

SPECIFICATIONS:

BLOWER DIMS:	13"D x 13"W x 13"H
OVERALL WEIGHT:	Approx. 30 lbs.
BLOWER MATERIAL:	16 ga. Carbon Steel w/ Epoxy Coating on All Surfaces
ELECTRICAL:	115/1/60, Approximately 2 Amps
SOUND LEVEL:	58 dba to 66 dba (Low to High)
CFM:	Pre-filter & 99.97% HEPA Filter: 350 Pre-filter & Heavy Duty Carbon(10 lbs.): 300
WARRANTY:	Limited one year warranty from date of shipment on defects due to materials or workmanship.
PATENTED DESIGN:	USA PATENT#: 5,843,197

All Specifications are nominal and subject to change without notice.
Not intended explosive or highly combustible applications.